

# Enhancing surgical outcomes: the vital role of continuous support from surgical material suppliers

Bruno Campos de França<sup>1</sup>, Artur Henrique Galvão Bruno da Cunha<sup>2</sup>, Fábio Coelho Alves Silveira<sup>2</sup>,  
Marcelo Moraes Valença<sup>2</sup>

<sup>1</sup>MedCirúrgica, Recife, Pernambuco, Brazil

<sup>2</sup>Unimed Recife, Pernambuco, Brazil



mmvalenca@yahoo.com.br

Edited by:  
Fernando Cruz

The advancement of neurosurgery in recent decades has been propelled by the surgeon's ability to recognize challenges during access to anatomical structures or tumor removal and to think creatively, leading to innovation in the production of instruments that facilitate surgical procedures (1-3). This innovation aims to increase surgical speed, reduce the risk of inadvertent patient injuries, and enable the surgical removal of tumors or treatment of conditions without increasing the chance of permanent functional damage. With the progress of advanced technological possibilities and greater interdisciplinary exchange, many new surgical instruments are entering the market, and the hands of surgeons are eager for new options that facilitate their patient care duties.

One of the surgical procedures that has greatly benefited from technological advances is endoscopic surgery for intra-sellar tumors and intraventricular neuroendoscopic surgery. Surgeons typically use an endonasal transsphenoidal approach to reach the base of the sella turcica, located just a few millimeters between the internal carotid arteries. Or coagulation of the choroidal plexus during treatment of hydrocephalus. For instance, specialized surgical instruments like the PB SafeBlend 450 are crucial for accessing structures deep within the body, several centimeters away from the body surface. With its bipolar coagulation mechanism, it can effectively coagulate bleeding foci without causing damage to neighboring structures, resulting in minimal thermal damage. There are two models, one for transsphenoidal use and another bayonet type (angled shaft) for microsurgery use with a microscope (Figures 1 and 2). Note that the tip has a grasping movement and is silver-coated to prevent adherence to the coagulation tissue.

Submitted: March 15, 2024  
Accepted: March 30, 2024




Figure 1. Illustration of two models PB SafeBlend - one designed for transsohenoidal (1) use and another featuring a bayonet-style (angled shaft) for microsurgical procedures under a microscope (2). The tip exhibits a grasping movement and is silver-coated to prevent tissue adherence during coagulation.

Therefore, we conclude that ongoing support from surgical material suppliers is fundamentally important to ensure safe surgery, with a lower chance of complications, fewer hospitalization days, and a minimal reoperation rate. This leads to lower mortality rates and higher patient satisfaction, ultimately resulting in lower costs for the patient, hospital, and health insurance plan in the long term.

Artur Henrique Galvão Bruno da Cunha  
<https://orcid.org/0000-0002-8589-8563>  
 Marcelo Moraes Valença  
<https://orcid.org/0000-0003-0678-3782>



Figure 2. Illustration depicting the use of instrumentation during an intraventricular procedure within the cerebrospinal fluid in a patient with hydrocephalus, aimed at coagulating the choroid plexus (observe the use of the forceps PB SafeBlend 450 to coagulate the vessel between the pinching movement). The tip exhibits a grasping movement and is silver-coated to prevent tissue adherence during coagulation.

## References

1. Weinreich HM, Jin V, Crowell L, Skovlund SM, Williams QL, Buhimschi IA. Surgical Instrument Designers and Inventors-Where are the Women? *Am Surg.* 2023 Nov;89(11):4835-41.
2. Bohl MA, Baranoski JF, Sexton D, Nakaji P, Snyder LA, Kakarla UK, et al. Barrow Innovation Center Case Series: Early Clinical Experience with Novel Surgical Instrument Used To Prevent Intraoperative Spinal Cord Injuries. *World Neurosurg.* 2018 Dec;120:e573-e9.
3. White EJ, Cunnane EM, McMahon M, Walsh MT, Coffey JC, O'Sullivan L. Mechanical characterisation of porcine non-intestinal colorectal tissues for innovation in surgical instrument design. *Proc Inst Mech Eng H.* 2018 Aug;232(8):796-806.

**Conflict of interest:** This article was commissioned by Medcirúrgica, a Razek partner, and the national distribution of the product is carried out by Medcirúrgica, Brasil.